

TRUPOWER™ Enzymatic Manual Detergent & Soak

A 2-in-1 product that makes manual warewashing fast, easy
and cost effective



Easy-to-use:

A multi-purpose pot & pan detergent plus soak simplifies procedures, product selection and training



Cleans effortlessly:

Enzymatic solution eliminates hard-to-remove soils and makes manual warewashing easy



Controls cost:

Ready-to-use dispensing helps prevent overuse



Reduces labor & waste:

Ultra-concentrated formula minimizes plastic packaging & results in fewer change outs

NEW



50%
less scrubbing¹

¹ Based on Ecolab internal testing compared to Pantastic Pot & Pan Detergent.

Cost saving, easy to use and reduces waste



Cost-efficient

Ultra concentrated solution; closed-loop dispensing prevents overuse compared to wasteful free pouring.




Easy to use

Multi-purpose, enzyme enhanced, detergent & soak built to tackle soils from soak to scrub to shine.

Reduces
plastic
waste by **96%**


when switching from Ecolab ready-to-use detergent¹

**One
2-gal bag**



=

**38
bottles**



of TRUPOWER Enzymatic Manual Detergent & Soak

of 32 oz RTU detergent¹



Ordering information

PRODUCT	SKU	PACK SIZE
TRUPOWER™ Enzymatic Manual Detergent & Soak	6102600	1 - 2 gal

¹ Based on Ecolab's dilution rate for TRUPOWER Enzymatic Manual Detergent & Soak compared to Ecolab's ready-to-use Pantastic Detergent bottle.

Contact your Ecolab Representative or call
1 800 35 CLEAN for more information

